

DIGITAL HAND-HELD "POCKET" REFRACTOMETER

65th
anniversary

PAL-α

Cat.No.3840

Special Anniversary model
with a wide Brix range **0 to 85%!!**



Actual
Size

Measurement Method



Apply 2-3 drops onto the prism surface.



Press the START key.



The measurement value will be displayed in 3 seconds.

Application

- * Outdoors in the vegetable fields and fruit orchards
- * Indoors for beverages (fruits juice, coffee, soft drink, etc) and foods (soy sauce, sauce, seasoning, soup, jam, marmalade, maple syrup, etc.)
- * At fruit and vegetable packing facilities and in the supermarket
- * At restaurants or shops
- * At junior high schools, high schools, or universities as course material

Features

- * Calibration with water only
- * Measurement in 3 seconds
- * Extremely water resistant (IP65)
- * Light & Compact, 100g
- * Ergonomically designed for easy, one-handed operation
- * Convenient storage case
- * Automatic Temperature Compensation (ATC)
- * Designed with a strap holder option



Specifications

Measurement range	: Brix 0.0 to 85.0%
Resolution	: Brix 0.1%
Measurement accuracy	: Brix $\pm 0.2\%$
Repeatability	: Brix $\pm 0.1\%$
Ambient Temperature	: 10 to 40°C
Measurement Temperature	: 10 to 75°C (Automatic Temperature Compensation)
Sample volume	: 0.3ml or more
Measurement time	: 3 seconds
Power supply	: AAA alkaline battery $\times 2$
Battery life	: Approx. 11,000 measurements (Using new alkaline batteries)
International Protection Class	: IP65 Water resistant
Dimensions and weight	: 55(W) $\times$ 31(D) $\times$ 109(H)mm, 100g

All ATAGO refractometers are designed and manufactured in Japan.

HACCP GMP GLP

ATAGO products comply with HACCP, GMP, and GLP system standards.

ATAGO CO., LTD.

Headquarters: The Front Tower Shiba Koen, 23rd Floor
2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan
TEL : 81-3-3431-1943 FAX : 81-3-3431-1945
overseas@atago.net http://www.atago.net/

ATAGO U.S.A., Inc.

TEL : 1-425-637-2107 customerservice@atago-usa.com

ATAGO INDIA Instruments Pvt. Ltd.

TEL : 91-22-2833-8038 / 8076 customerservice@atago-india.com

ATAGO (THAILAND) Co., Ltd.

TEL : 662-982-8718-9 customerservice@atago-thailand.com

ATAGO BRASIL Ltda.

TEL : 55 16 3916-6000 customerservice@atago-brasil.com

ATAGO ITALIA s.r.l.

TEL : 39 2 36557267 customerservice@atago-italia.com

ATAGO CHINA Guangzhou Co., Ltd.

TEL : 86 20-38106065 info@atago-china.com



* Specifications and appearance are subject to change without notice.

V.02 1207K Printed in Japan